Specifications

Low Flux Hollow Fiber Haemodialyser

Type	UF Coefficient (ultrafiltration) (ml/mmHg.h)	KOA ml/min	Clearance (ml/min)										Blood priming volume (ml)	Surface Area (m ²)
			Urea		Creatinine		Phosphate		VitaminB ₁₂		Inulin			
			Q _B =200	Q _B =300	Q _B =200	Q _B =300	Q _B =200	Q _B =300	Q _B =200	Q _B =300	Q _B =200	Q _B =300		
PSR 10 LF	7.3	698	187	238	167	201	143	161	79	84	-	-	59	1.0
PSR 13 LF	9.5	824	195	250	180	223	153	182	93	100	-	-	69	1.3
PSR 14 LF	11.9	977	196	253	184	226	157	188	96	117	-	-	74	1.4
PSR 16 LF	13.9	1240	197	274	189	243	164	196	120	135	-	-	86	1.6
PSR 18 LF	16.7	1570	198	284	189	265	174	215	139	156	-	-	105	1.8
PSR 20 LF	17.8	1610	199	285	191	268	177	220	145	161	-	-	107	2.0
PESR 10 LF	9.3	767	186	245	182	221	158	182	104	113	-	-	59	1.0
PESR 13 LF	12.2	961	190	260	184	235	163	205	131	139	-	-	72	1.3
PESR 16 LF	14.8	1292	194	276	186	248	170	219	142	153	-	-	89	1.6
PESR 18 LF	22.6	1450	198	281	189	262	174	226	147	161	-	-	110	1.8
PESR 20 LF	25.5	1662	199	286	196	270	186	242	152	174	-	-	114	2.0

High Flux Hollow Fiber Haemodialyser

Type	UF Coefficient (ultrafiltration) (ml/mmHg.h)	KOA ml/min	Q _B =200	Q _B =300	Q _B =200	Q _B =300	Q _B =200	Q _B =300	Q _B =200	Q _B =300	Q _B =200	Q _B =300	Blood priming volume (ml)	Surface Area (m ²)
PSR 100 HF	34.6	930	190	258	173	217	160	199	113	127	78	84	59	1.0
PSR 130 HF	46.4	1030	195	264	180	234	174	219	130	146	93	103	69	1.3
PSR 140 HF	50.6	1190	196	272	186	247	179	230	139	159	104	113	74	1.4
PSR 160 HF	59.4	1570	197	284	193	267	188	249	153	178	117	130	86	1.6
PSR 180 HF	63.7	1660	198	286	194	276	192	259	161	190	125	141	105	1.8
PSR 200 HF	74.6	1900	199	290	196	283	195	277	168	208	131	152	107	2.0
PESR 130 HF	58.3	1140	196	270	192	255	181	236	147	168	107	116	72	1.3
PESR 160 HF	67.0	1610	197	285	193	269	189	254	160	186	120	134	89	1.6
PESR 180 HF	84.2	1770	198	288	197	279	194	262	165	194	127	144	110	1.8
PESR 200 HF	87.5	1980	199	291	197	286	196	279	172	211	132	154	114	2.0

Specifications and performance data at QB: 200/300 ml/min, QD: 500 ml/min, QF: 0 ml/min, T: 37°C
Performance data were measured in vitro according to standard BS ISO 8637-1:2024

Typical values obtained with an individual batch of fiber, the clinical use may determine difference in results in relation to different Ultrafiltration and measuring technique possible verification between batches of fibers.

UF measurement using Bovine/Human blood (Hct 32%, Protein 60 g/l)

The tolerance for clearance numbers is $\pm 10\%$

HF: High Flux LF: Low Flux

Direction for use

Intended Purpose: Hollow-fiber haemodialysis dialyser is a filter that functions as an artificial kidney and typically used in a haemodialysis system to remove impurities/fluid from the blood of a patient. It typically consists of cylindrical containers with thousands of longitudinally arranged hollow-fiber capillary tubes (e.g., polymer, modified cellulose) through which the blood flows. The tube walls function as a semi-permeable membrane, permitting passage of larger molecules from the blood to the dialysate on the outside of the tubes for removal. This is a single-use device.

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**contraindications**

In general, dialysis disease is not sex-dependent. PS and PES haemodialysers of meditechsys are not recommended for pediatric, newborn and toddler use. Factors such as diabetes, hypertension, old age and etc. Increase the likelihood of this disease. Dialyser choice is according to its surface area, weight of patient and patient clinical test (toxins concentration) by physician.

There are no specific contraindications for the dialyser itself. Avoid use in patients with known hypersensitivity to particular dialyser (PS and PES). Usage of haemodialyser is contraindicated in patients presenting with a known allergy to heparin or who have type II thrombocytopenia caused by heparin (HIT Syndrome type II).

PS (Polysulfone series) and PES (Polyethersulfone series) haemodialyzers are non-pyrogenic devices and sterilized with gamma irradiation.

Warning

Do not use for any other purposes than haemodialysis.

- The dialyser should only be used as directed by a physician.
- If foaming, blood leakage, blood coagulation and hemolysis occur during the use of this product, take appropriate measures according to a physician's instructions.
- DO NOT REUSE THIS DEVICE.

If reuse happened, risk clotting, severe infection and diminished performance is expected.

- If the patient exhibits any abnormal symptoms such as discomfort, pruritus, urticaria, asthmatic reaction, hypertension and/or arrhythmia during the use of this product, take appropriate measures according to a physician's instructions.
- Dialyser contact with chemicals and disinfections affect the properties of housing, polyurethane and fiber. The operating safety does not longer exist and the manufacturer's liability is excluded.
- Heparinize the extracorporeal circuit during the priming, these dosage of heparin and its mode of administration during dialysis are established by attending physician. The coagulation time should be checked regularly.

Device performance is diminished if it is used below mentioned flow rates, below mentioned pressure or in particular orientations. This instruction should be followed.

NB: Meditechsys will not be responsible for injury of patient or any person or damage to any object that is attributed to bad transport, bad storage, shipping accidents, abnormal use, user error and reasonably foreseeable misuse.

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**Caution**

It is recommended to use of a dialysis machine equipped with a volumetric ultrafiltration unit, an accurate UF control system and a precise weigh system for haemodialysis.

1. Caution before use

- Do not use the device if package is damaged or if protective end caps are not in place.
- Avoid air intake and contamination during rinsing and priming operations.
- Blood-side: Rinse extracorporeal circuit with sterile physiological saline (recommended 1000 ml and 500 ml of heparinized physiological saline, 2000 U/500ml) to minimize the side effect of sterilization residues or to prepare the fibers for blood flow at a flow rate of 100 ml/min.
- Perform a leak test to check the integrity of the bloodline and dialyser.
- Start dialysis mode immediately after priming operation.

2. Caution in use

- Continuously monitor the arterial and venous pressures in the bloodline and check for blood leakage during dialysis.
- Set TMP alarm (max. 500 mmHg).

3. Caution after use

- Dispose of the dialyser immediately after use by suitable means as contaminated medical waste and according to any prevailing environmental regulations.
- The weight loss should be monitored.

4. Caution for storage

- Store at 5°C~35°C avoiding exposure to direct sunlight, severe vibration and relative humidity between %25 - %65.

5. Connectors

- The blood ports of haemodialysers comply with the requirements of BS ISO 8637-1: 2024. They can connect to female lure lock connectors of blood line set.
- The dialysis fluid port of haemodialyser comply with the requirements of BS ISO 8637-1: 2024. They can connect to the dialysis fluid port of standard dialysis machines.

Rinsing and air removal

1. Unpack the dialyser and clip onto the holder ensuring that the venous side is directed upward and the arterial side is directed downward.
2. Connect the arterial line to the physiological

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

3. Connect aseptically the arterial and venous line to the dialyser and remove air completely by running physiological saline (recommended volume 1000ml), at a flow rate of about 100ml/min.

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4. After confirming that there are no bubbles in arterial and venous headers, stop the blood pump, turn the dialyser 180° to allow removal of bubbles from dialysate. Connect the dialysate lines to the dialyser and run dialysate at a flow of about 500ml/min under dialysate pressure of 0 mmHg through the extracorporeal circuit for about 5 minutes.

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
Priming

Following the rinsing and air removal, start the blood pump, run about 500 ml of physiological saline (eventually heparinized, 2000 IU/500 ml) through the arterial line at a flow rate of about 100 ml/min; Prime the dialyser with heparinized saline and dialysate flow.

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
Leakage test

- Activate the bypass valve on the machine and connect the dialyser lines to the machine.
- Prime the arterial, venous line and the dialyser completely with physiological saline by operating the blood pump, then stop the pump operation.
- Clamp the arterial line near the dialyser and the distal end of the venous line with forceps.
- Place the clamped distal end about 1m below the dialyser and remove the forceps. (This results in application of a negative pressure of about 70 mmHg to the blood compartment of the dialyser.) Examine for any or no continuous stream of bubbles in the venous header to check for leakage from the dialyser; if observed, replace the dialyser with a new one.

Start of dialysis

1. Stop the blood flow. Prepare the blood access site and connect to the arterial line. Remove the forceps from the arterial and venous lines. While running dialysate at a flow of about 500 ml/min, operate the blood pump at a flow rate of about 100 ml/min.
2. Confirm that no air bubbles remain in venous header or venous bloodline.


(Fig.4)

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3. Fully prime the arterial and venous lines including the dialyser with blood by operating the blood pump, and then stop the pump operation. Clamp the distal end of the venous line with forceps.
4. Prepare the blood return site and connect to the venous line. After confirming that there are no bubbles in the line, remove the forceps from the line. After checking that there are no forceps on the lines and no line folding, operate the blood pump at a low flow rate. Avoid excessive pressure to the lines and the dialyser to prevent leakage from the dialyser and separation of each of the connection.
5. Increase the blood pump rate to 200 or 300 ml/min according to patient condition.

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**Operations during dialysis**

1. If stopping the blood pump is required during dialysis due to insufficient blood flow or other reasons, lower the Ultrafiltration rate UF to a safe minimum.

2. Set UF rate carefully to avoid excessive water removal. Reduce the blood flow rate accordingly to avoid a disequilibrium syndrome.

3. In case of blood leakage alarm, stop the dialysate supply and recover blood, then replace the dialyser with a new one.

Dialysis treatment-end and blood recovery

1. Lower the UF and TMP as low as possible. Stop the blood pump, clamp the arterial line and remove the line from the arterial blood access site; then connect the line to the physiological saline vial for blood recovery.

2. Unclamp the arterial line and run 100~200ml of physiological saline at a flow rate of about 100 ml/min to expel blood from the arterial and venous lines and the dialyser.

3. After blood recovery, discard the blood line sets and the dialyser by suitable means as medical waste.

Caution:

In Long-term treatment, chronic side effects such as amyloidosis, anemia, depression, thrombocytopenia, coronary arteries, and congestive heart diseases may occur, you must consult your doctor.

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Read the label on the product or on the carton prior to use:

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Typical values obtained with an individual batch of fiber, the clinical use may determine difference
in results in relation to different Ultrafiltration and measuring technique possible verification
between batches of fibers.

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ICD GROUP

Meditechsyst

Manufacture Address: Meditechsyst Co.

Head Office: South Apt. 2, 5th floor, No. 44, East 15 Mirzaei Alley, Saadatabad St., Tehran, IRAN

Postal Code: 1998898639

Tel: +98(21) 22074001 **Fax:** +98(21) 22073101

Factory: Plots 295 and 296, 3rd Tangestan St., West Hafez Blvd., Aboureyhan Blvd., Eshtehard Industrial Park, Eshtehard County, Alborz Province, IRAN

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Email: eurep@nephrocan.com

Note: Meditechsyst is not responsible for any shipping accidents, bad storage and abnormal use.

Any serious incident that has occurred in relation to the device should be reported immediately to the manufacturer (Meditechsyst Co.), its European Authorized Representative, and the competent authority of the Member State in which the device has been made available (imported) or the incident occurred, and/or where the user and/or patient is established.

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NephroCan GmbH. & Co. KG