Specifications

Low Flux Hollow Fiber Haemodialyser														
Type	UF Coefficient (ultrafiltration) (ml/mmHg.h)	KOA ml/min	Clearance (ml/min)										Blood priming volume (ml)	Surface Area (m ²)
			Urea		Creatinine		Phosphate		VitaminB ₁₂		Inulin			
			Q _B = 200	Q _B = 300	Q _B = 200	Q _B = 300	Q _B = 200	Q _B = 300	Q _B = 200	Q _B = 300	Q _B = 200	Q _B = 300		
PS 10 LF	6.8	637	183	231	164	196	140	158	73	78	-	-	59	1.0
PS 13 LF	8.8	746	191	243	176	218	150	178	86	93	-	-	69	1.3
PS 14 LF	10.6	861	193	253	179	226	158	184	95	108	-	-	74	1.4
PS 16 LF	12.9	1064	195	266	184	237	161	192	111	125	-	-	86	1.6
PS 18 LF	15.5	1292	197	276	188	259	171	211	129	144	-	-	105	1.8
PS 20 LF	16.5	1320	198	277	190	261	174	216	134	149	-	-	107	2.0
PES 10 LF	8.9	698	184	238	180	215	154	177	96	105	-	-	59	1.0
PES 13 LF	11.6	861	188	253	182	229	159	199	121	129	-	-	72	1.3
PES 16 LF	14.1	1103	192	268	184	242	166	213	132	142	-	-	89	1.6
PES 18 LF	21.5	1214	197	273	187	255	170	220	136	149	-	-	110	1.8
PES 20 LF	24.3	1351	198	278	194	263	181	235	141	161	-	-	114	2.0
High Flux Hollow Fiber Haemodialyser														
PS 100 HF	32	778	184	246	168	205	156	186	105	118	72	78	59	1.0
PS 130 HF	43	836	189	251	175	221	170	205	120	135	86	95	69	1.3
PS 140 HF	47	961	192	260	181	233	176	215	131	148	97	105	74	1.4
PS 160 HF	55	1145	195	270	191	252	183	233	142	165	108	120	86	1.6
PS 180 HF	59	1265	196	275	193	260	187	242	149	176	116	131	105	1.8
PS 200 HF	69	1380	198	279	195	267	190	259	156	193	121	141	107	2.0
PES 130 HF	54	916	191	257	186	241	177	221	136	156	99	107	72	1.3
PES 160 HF	62	1167	194	271	191	254	184	237	148	172	111	124	89	1.6
PES 180 HF	78	1321	198	277	195	263	189	245	153	180	118	133	110	1.8
PES 200 HF	81	1420	199	280	196	270	192	261	159	195	122	143	114	2.0

Specifications and performance data at QB: 200 or 300 ml/min, QD: 500 ml/min, QF: 0 ml/min, T: 37°C
Performance data were measured in vitro according to standards BS ISO 8637-1:2024

Typical values obtained with an individual batch of fiber, the clinical use may determine difference in results in relation to different Ultrafiltration and measuring technique possible verification between batches of fibers.

UF measurement using Bovine/Human blood (Hct 32%, Protein 60 g/l)

The tolerance for clearance numbers is $\pm 10\%$

HF: High Flux LF: Low Flux

Direction for use

Intended Purpose: Hollow-fiber haemodialysis dialyser is a filter that functions as an artificial kidney and typically used in a haemodialysis system to remove impurities/fluid from the blood of a patient. It typically consists of cylindrical containers with thousands of longitudinally arranged hollow-fiber capillary tubes (e.g., polymer, modified cellulose) through which the blood flows. The tube walls function as a semi-permeable membrane, permitting passage of larger molecules from the blood to the dialysate on the outside of the tubes for removal. This is a single-use device.

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Contraindications: In general, dialysis disease is not sex-dependent.

PS and PES haemodialysers of meditechsys are not recommended for pediatric, newborn and toddler use. Factors such as diabetes, hypertension, old age and etc. Increase the likelihood of this disease. Dialyser choice is according to its surface area, weight of patient and patient clinical test (toxins concentration) by physician.

There are no specific contraindications for the dialyser itself. Avoid use in patients with known hypersensitivity to particular dialyser (PS and PES) or Ethylene Oxide (EO) gas. Usage of haemodialyser is contraindicated in patients presenting with a known allergy to heparin or who have type II thrombocytopenia caused by heparin (HIT Syndrome type II).

Warning

- Do not use for any other purposes than haemodialysis.

- The dialyser should only be used as directed by a physician.

- If foaming, blood leakage, blood coagulation and hemolysis occur during the use of this product, take appropriate measures according to a physician's instructions.

- DO NOT REUSE THIS DEVICE.

If reuse happened, risk clotting, severe infection and diminished performance is expected.

- If the patient exhibits any abnormal symptoms such as discomfort, pruritus, urticaria, asthmatic reaction, hypertension and/or arrhythmia during the use of this product, take appropriate measures according to a physician's instructions.

- Dialyser contact with chemicals and disinfections affect the properties of housing, polyurethane and fiber. The operating safety does not longer exist and the manufacturer's liability is excluded.

- Heparinize the extracorporeal circuit during the priming, these dosage of heparin and its mode of administration during dialysis are established by attending physician. The coagulation time should be checked regularly.

Device performance is diminished if it is used below mentioned flow rates, below mentioned pressure or in particular orientations. This instruction should be followed.

NB: Meditechsys will not be responsible for injury of patient or any person or damage to any object that is attributed to bad transport, bad storage, shipping accidents, abnormal use, user error and reasonably foreseeable misuse.

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Caution

It is recommended to use of a dialysis machine equipped with a volumetric ultrafiltration unit, an accurate UF control system and a precise weigh system for haemodialysis.

1. Caution before use

- Do not use the device if package is damaged or if protective end caps are not in place.
- Avoid air intake and contamination during rinsing and priming operations.
- Blood-side: Rinse extracorporeal circuit with sterile physiological saline (recommended 1000 ml and 500 ml of heparinized physiological saline, 2000 U/500ml) to minimize the side effect of sterilization residues or to prepare the fibers for blood flow at a flow rate of 100 ml/min.
- Perform a leak test to check the integrity of the bloodline and dialyser.
- Start dialysis mode immediately after priming operation.

2. Caution in use

- Continuously monitor the arterial and venous pressures in the bloodline and check for blood leakage during dialysis.
- Set TMP alarm (max. 500 mmHg).

3. Caution after use

- Dispose of the dialyser immediately after use by suitable means as contaminated medical waste and according to any prevailing environmental regulations.
- The weight loss should be monitored.

4. Caution for storage

- Store at 5°C~35°C avoiding exposure to direct sunlight, severe vibration and relative humidity between %25 %65.

5. Connectors

The blood ports of haemodialysers comply with the requirements of BS ISO 8637-1: 2024. They can connect to female luer lock connectors of blood line set. The dialysis fluid port of haemodialyser comply with the requirements of BS ISO 8637-1: 2024. They can connect to the dialysis fluid port of standard dialysis machines.

Rinsing and air removal

- Unpack the dialyser and clip onto the holder ensuring that the venous side is directed upward and the arterial side is directed downward.
- Connect the arterial line to the physiological saline vial.
- Connect aseptically the arterial and venous line to the dialyser and remove air completely by running physiological saline (recommended volume 1000ml), at a flow rate of about 100ml/min.

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

4. After confirming that there are no bubbles in arterial and venous headers, stop the blood pump, turn the dialyser 180° to allow removal of bubbles from dialysate. Connect the dialysate lines to the dialyser and run dialysate at a flow of about 500ml/min under dialysate pressure of 0 mmHg through the extracorporeal circuit for about 5 minutes.

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
Priming

Following the rinsing and air removal, start the blood pump, run about 500 ml of physiological saline (eventually heparinized, 2000 IU/500 ml) through the arterial line at a flow rate of about 100 ml/min; Prime the dialyser with heparinized saline and dialysate flow.

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Leakage test

- Activate the bypass valve on the machine and connect the dialyser lines to the machine.
- Prime the arterial, venous line and the dialyser completely with physiological saline by operating the blood pump, then stop the pump operation.
- Clamp the arterial line near the dialyser and the distal end of the venous line with forceps.
- Place the clamped distal end about 1m below the dialyser and remove the forceps. (This results in application of a negative pressure of about 70 mmHg to the blood compartment of the dialyser.) Examine for any or no continuous stream of bubbles in the venous header to check for leakage from the dialyser; if observed, replace the dialyser with a new one.

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
Start of dialysis

1. Stop the blood flow. Prepare the blood access site and connect to the arterial line. Remove the forceps from the arterial and venous lines. While running dialysate at a flow of about 500 ml/min, operate the blood pump at a flow rate of about 100 ml/min.
2. Confirm that no air bubbles remain in venous header or venous bloodline.

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
3. Fully prime the arterial and venous lines including the dialyser with blood by operating the blood pump, and then stop the pump operation. Clamp the distal end of the venous line with forceps.
4. Prepare the blood return site and connect to the venous line. After confirming that there are no bubbles in the line, remove the forceps from the line. After checking that there are no forceps on the lines and no line folding, operate the blood pump at a low flow rate. Avoid excessive pressure to the lines and the dialyser to prevent leakage from the dialyser and separation of each of the connection.
5. Increase the blood pump rate to 200 or 300 ml/min according to patient condition.

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Operations during dialysis

1. If stopping the blood pump is required during dialysis due to insufficient blood flow or other reasons, lower the Ultrafiltration rate UF to a safe minimum.
2. Set UF rate carefully to avoid excessive water removal. Reduce the blood flow rate accordingly to avoid a disequilibrium syndrome.
3. In case of blood leakage alarm, stop the dialysate supply and recover blood, then replace the dialyser with a new one.

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Dialysis treatment-end and blood recovery

1. Lower the UF and TMP as low as possible. Stop the blood pump, clamp the arterial line and remove the line

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Meditechsys

Instruction for Use of Hollow Fiber Haemodialyser

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JAM TRITA
GROUP

ICD GROUP

from the arterial blood access site; then connect the line to the physiological saline vial for blood recovery.

2. Unclamp the arterial line and run 100~200ml of physiological saline at a flow rate of about 100 ml/min to expel blood from the arterial and venous lines and the dialyser.

3. After blood recovery, discard the blood line sets and the dialyser by suitable means as medical waste.

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Caution:

In Long-term treatment, chronic side effects such as amyloidosis, anemia, depression, thrombocytopenia, coronary arteries, and congestive heart diseases may occur, you must consult your doctor.

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Read the label on the product or on the carton prior to use:

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Meditechsys

Manufacture Address: Meditechsys Co.

Head Office: South Apt. 2 ,5th floor, No. 44, East 15 Mirzaei Alley, Saadatabad St., Tehran, IRAN

Postal Code: 1998898639

Tel: +98(21) 22074001 **Fax:** +98(21) 22073101

Factory: Plots 295 and 296 ,3rd Tangestan St., West Hafez Blvd., Aboureyhan Blvd., Eshtehard Industrial Park, Eshtehard County, Alborz Province, IRAN

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Email: eurep@nephrocan.com

Note: Meditechsys is not responsible for any shipping accidents, bad storage and abnormal use.

Any serious incident that has occurred in relation to the device should be reported immediately to the manufacturer (Meditechsys Co.), its European Authorized Representative, and the competent authority of the Member State in which the device has been made available (imported) or the incident occurred, and/or where the user and/or patient is established.

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NephroCan GmbH. & Co. KG