1. Intended Purpose

The Bloodline set and Fistula Needle are a sterile, single-use medical device intended to provide a conduit for the extracorporeal circulation of blood between a patient's vascular access and a hemodialyzer during hemodialysis treatment. The device is intended for use exclusively by trained healthcare professionals in clinical or dialysis center settings.

2. Indication

Used in the treatment of patients with acute or chronic renal failure, when hemodialysis is prescribed by a physician. The Bloodline Set and Fistula Needle are non-pyrogenic, phthalate-free, and sterilized either by Ethylene Oxide (EO) gas or gamma irradiation. Please refer to the product packaging for information regarding the sterilization method.

3. Contraindications

Patients with hypersensitivity to Ethylene Oxide gas (EO) should not use products sterilized with this gas.

4. Side effect

The Bloodline set and the fistula needle themselves have no pharmacological effects; any reported side effects are generally related to their clinical use or specific patient reactions.

5. Warning and Cautions

1. Do not use the product if the packaging or device is damaged, or if caps are missing.
2. When using the Bloodline set for HDF, high-flux filters must be used.
3. Ensure the fistula needle cap is securely in place.
4. The product is intended only for hemodialysis and must not be used for other applications.
5. Do not reuse or re-sterilize the product.
6. Kinking of the tubing may cause hemolysis; discontinue use if observed.
7. Ensure all connections are properly secured to prevent blood leakage or air entry.
8. Do not use air flush methods to return blood to the patient at the end of dialysis.
9. After use, dispose of the product immediately as medical waste in accordance with local regulations.
10. Only trained healthcare personnel are authorized to use this product.
11. The expiration date applies to a product that has been properly stored and remains intact.
12. Pharmed Medical Industries is not responsible for improper storage or transport, misuse, or use of extracorporeal blood circuits with incompatible equipment.
13. Avoid excessive pressure or abnormal blood flow, as this may lead to hemolysis or device damage.
14. Ensure all air is completely removed from the system before starting dialysis to prevent air embolism.
15. Do not connect the Bloodline set to incompatible or non-approved dialysis machines.
16. The pediatric bloodline set is recommended for children weighing over 12 kg.

6. Direction for use

1. Check the integrity of the packaging and confirm that the product is sterile.
2. Remove the Bloodline set from its packaging.
3. Use the product immediately after opening.
4. Connect the arterial line of the set to a normal saline solution.
5. If the set includes a drainage bag at the end of the venous line, hang it on the designated stand and ensure the clamp is open.

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6. Connect the arterial and venous lines of the set to the dialysis machine. Prime the lines completely by running normal saline through the set at approximately 100 ml/min (recommended priming volume: 1000 ml).

7. Priming:

After washing and priming, pass approximately 500 ml of normal saline through the arterial line at a flow rate of 100 ml/min. It is recommended to heparinize the Bloodline set during preparation. The heparin dosage and administration method for the dialysis session should be determined by the treating physician. Clamp the necessary parts of the arterial and venous lines. Complete all connections required to start dialysis. Check for any leaks (continuous bubbles) in the Bloodline set and dialyzer; if a leak is observed, replace the circuit. Figure 1. Also, for priming of the arterial and venous lines, refer to the instructions provided by the dialysis machine manufacturer.

8. Starting Hemodialysis Treatment

a) Prepare and disinfect the vascular access site using gloves. Insert arterial and venous needles at the appropriate sites with sufficient distance from the vessel wall. Open the clamp and cap of the arterial needle; once blood is observed in the tubing, reclamp it. Repeat the same procedure for the venous needle.

b) Connect the arterial line to the patient's vascular access. Remove clamps from arterial and venous lines. Set the dialysis solution flow rate to 500 ml/min and start the blood pump at 100 ml/min.

c) Ensure no air bubbles are present in the venous line.

d) Once the arterial, venous lines, and dialyzer are fully primed with blood, stop the pump. Clamp the end of the venous line and the drainage bag, and disconnect the venous line from the bag. Dispose of the bag according to the proper procedure.

e) Connect the venous line to the patient's bloodstream. After ensuring no air remains in the line, remove the clamp. Check that no other clamps are closed and that the lines are free of kinks. Start the blood pump gradually. Avoid applying excessive pressure on the lines or dialyzer to prevent leaks or disconnections.

f) After confirming no air bubbles in the arterial and venous lines, rotate the dialyzer 180° to remove any bubbles from the dialysis solution. If bubbles are still present, repeat the rotation. Set the blood pump to the prescribed speed and, with the venous side of the dialyzer elevated, run the flow for approximately 10 minutes. Note: This procedure applies to dialysis lines with a drain bag. For Bloodline sets without a drainage bag, after stopping the pump, simply clamp the end of the venous line.

9. End of Dialysis and Blood Recovery

a) Stop the blood pump. Clamp the arterial line and disconnect it from the vascular access, then connect it to normal saline.

b) Open the clamp on the arterial line and run approximately 100-200 ml of normal saline at 100 ml/min to wash out any remaining blood from the arterial line, venous line, and dialyzer.

▲ The maximum allowable transmembrane pressure (TMP) for the Bloodline set is 500 mmHg.

▲ Fistula needles of gauges 15, 16 and 17 are recommended for use in hemodialysis.

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
17G

Recommended Fistula Needle


Figure 1: Typical Fluid circuit diagram

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7. Recommended suitable of Hemodialysis machine

- Blood Line Set, SU – Recommended for use with most dialysis machines.
- Blood Line Set, SF – Optimized for Fresenius machines, but also compatible with most dialysis machines.
- Blood Line Set, SB – Recommended for use with B. Braun machines.
- Blood Line Set, SG – Recommended for use with Gambro machines.
- Blood Line Set, HDF Online – Recommended for use with Dialog+ HDF Online systems.
- Blood Line Set, Pediatric – Recommended for pediatric hemodialysis treatment and compatible with most dialysis machines.

▲ The selection of a suitable bloodline set should always be made by the dialysis nurse or healthcare professional, based on the specific dialysis machine in use.

Bloodline Half-Sets

Venous Line, SF and Arterial Line, SF

→ These half-sets are intended for situations where replacement of only part of the bloodline (arterial or venous) is required during dialysis, without the need to replace the entire set.

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Note: Only products marked with an asterisk (*) have the CE mark.

Fistula Needle Specification			RenaNeedle®	
No.	Model Name	Wing Type/Color	Product Code (REF.)	UDI-DI
Sterile with EO				
1	Arterial Fistula Needle ✱ 15G (1.8*25) 150 mm	Rotating/Blue	4218366	06265915400770
2	Venous Fistula Needle ✱ 15G (1.8*25) 150 mm		4218364	06265915400749
3	Fistula Needle Set (Venous & Arterial) ✱ 15G (1.8*25) 150 mm		4218384	06265915401036
4	Arterial Fistula Needle ✱ 16G (1.6*25) 150 mm	Rotating/Green	4218340	06265915400367
5	Venous Fistula Needle ✱ 16G (1.6*25) 150 mm		4218342	06265915400626
6	Fistula Needle Set (Venous & Arterial) ✱ 16G (1.6*25) 150 mm		4218302	06265915401074
7	Arterial Fistula Needle ✱ 17G (1.4*25) 150 mm	Rotating/Orange	4218354	06265915400374
8	Venous Fistula Needle ✱ 17G (1.4*25) 150 mm		4218356	06265915400633
9	Fistula Needle Set (Venous & Arterial) ✱ 17G (1.4*25) 150 mm		4218359	06265915401142
Sterile with Gamma Irradiation				
1	Arterial Fistula Needle-R 15G (1.8*25)150mm	Rotating/Blue	4218301	06265915401098
2	Venous Fistula Needle-R 15G (1.8*25)150mm		4218303	06265915401081
3	Fistula Needle Set-R (Venous & Arterial) 15G (1.8*25)150mm		4218304	06265915401043
4	Arterial Fistula Needle-R 16G (1.6*25)150mm	Rotating/Green	4218307	06265915401111
5	Venous Fistula Needle-R 16G (1.6*25)150mm		4218309	06265915401104
6	Fistula Needle Set-R (Venous & Arterial) 16G (1.6*25)150mm		4218311	06265915401050
7	Arterial Fistula Needle-R 17G (1.4*25)150mm	Rotating/Orange	4218313	06265915401128
8	Venous Fistula Needle-R 17G (1.4*25)150mm		4218316	06265915401135
9	Fistula Needle Set-R (Venous & Arterial) 17G (1.4*25)150mm		4218318	06265915401067

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Note :Only products marked with an asterisk (*) have the CE mark.

9. Storage Condition

Store at 5°C to 40°C, protected from direct sunlight, severe vibration, and at a relative humidity not exceeding 65%. Protect from freezing.

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NephroCan GmbH. & Co.KG