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A. Intended use

Concentrated Hemodialysis Acidic Solution that diluted with dialysis water and bicarbonate concentrate, yields dialysis fluid for use in dialysis.

Dialysis fluid is used in the dialysis process to remove waste products and metabolites from the blood of patients with kidney problems who are unable to eliminate toxins and excess fluids on their body.

B. Side Effects

Use of Concentrated Hemodialysis Acidic Solution may lead to some contraindications these are as follows: Hypotonia, nausea, vomiting and spasms have been reported.

C. Contraindications

- That the solution is not to be used for intravenous infusion.
- Use of Concentrated Hemodialysis Acidic Solution may lead to some contraindications these are as followed:
 1. Acidic Hemodialysis Concentrate I don't use Hypokalemia (applicable for potassium-free concentrate), and unmanageable coagulation disorders.
 2. Acidic Hemodialysis Concentrate III don't use Hyperkalemia (applicable for potassium containing concentrates).
 3. For patients with haemodynamic instability, an alternative method of extracorporeal treatment may be indicated.

D. Caution and Warning

1. single use.
2. Concentrated Hemodialysis Acidic Solution must be diluted with Purified water immediately before use in accordance with the relevant national requirement. Discard residues.
3. Always check the concentration of the ready to use dialysis solution. Use only clear solution.
4. Do not use from damaged containers.
5. Concentrates are Non-Sterile and non-pyrogenic.
6. In case of treatment interruptions (Multiple Hemodialysis sessions); make sure that the container of the concentrated Hemodialysis solution is covered well, and keep it clean and away from any source of contamination.
7. For product performance; refer to the concentration of Hemodialysis solution mentioned on the product label.
8. To be used only by prescription.
9. Pharmed medical industries will not be responsible for bad storage, bad transportation, abnormal usage and use of the product with any incompatible device.
10. Discard the product, immediately after use by suitable means as contaminated medical waste and according to any prevailing environmental regulations.

E. Storage

Store Between (4°C-30°C) avoiding exposure to direct sunlight, severe vibration and keep out of freezing (Short-term exposure to 40°C will not harm to acid concentrated)

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H. Manufacture

Name: Pharmed Medical Industries

Factory: 2nd Bustan, West Hafez Blvd., Eshtehard Industrial Zone, IRAN Office: 2nd floor, No. 44, Saadat Abad Ave., 1998898638, Tehran, IRAN Customer Care Services: +989123346849 +98(21)22143509

Telefax: +98 (26) 3777 7620-5 Telefax: +98 (21) 2207 3102-3

Email: info@pharmed.co Website: www.pharmed.co